**CitraFlow™**  
 Prefilled Syringe with Anticoagulant and Antimicrobial  
 Sodium Citrate Solution 46.7% (w/v)

**INDIRECT IN-LINE ANTICOAGULANT INDICATIONS AND CLINICAL USE**

**CitraFlow™** is indicated for use in maintaining patency of Vascular Access Devices (VAD's).

**CONTRAINDICATIONS**  
 Patients with known allergies or hypersensitivity to constituents. For a listing see the PRODUCT DESCRIPTION AND PACKAGING section.

**WARNINGS AND PRECAUTIONS**  
 Use aseptic technique. Single use only. Do not reuse. Not for direct intravenous injection. Do not use if package is opened or damaged. Do not use if the cap on the syringe is not intact. Visually inspect the contents of each prefilled syringe for clarity, particulate matter, precipitate, discoloration and leakage prior to use. Do not use if any of the aforementioned is observed. Do not allow air to be trapped in fluid path. Do not resterilize. Sterile unless package opened or damaged. May be placed on a sterile field.

**ADVERSE REACTIONS**  
 No known adverse reactions have been reported when the product is used as indicated. Paraesthesia and/or dysgeusia may occur if the product unintentionally passes into the vein.

**INSTRUCTIONS FOR USE USE ASEPTIC TECHNIQUE.**

- Prior to initiation of a venous access session, aspirate and discard the solution that was instilled in the Vascular Access Device (VAD) during the previous session in accordance with institution's protocol.
- At the completion of each venous access session, flush the lumens of the VAD with sterile 0.9% Sodium Chloride solution to remove any blood, medication or other substance remaining in the VAD.
- Aseptically open package containing **CitraFlow™** prefilled syringe(s). Inspect the contents of the syringe(s) according to WARNINGS AND PRECAUTIONS.
- Remove and discard syringe cap. Expel air from syringe. Do not allow air to be trapped in fluid path. Connect the syringe to the VAD. Use only with compatible luer lock connectors.
- Slowly, over a 5 to 10 second period, instill the exact priming volume of **CitraFlow™** solution into each VAD lumen, as a locking agent, in volumes corresponding to the luminal capacity as specified by the VAD manufacturer.
- Remove and discard unused portions and empty syringes according to the institution's biohazardous waste policy.

**ACTION**  
 In the extracorporeal circuit, Sodium Citrate causes anticoagulation by chelating ionized Calcium into a soluble complex. Calcium is an integral ion involved in the clotting cascade. Local removal by Citrate prevents the activation of clotting cofactors, factor X and prothrombin and the ultimate formation of fibrin. When used as indicated, **CitraFlow™** eliminates the risks of systemic anticoagulation associated with the use of heparin locking solutions. The antimicrobial effect is through binding and removal of Ca2+ in the surrounding area. Ca2+ may regulate certain genes responsible for growth and survival of microbes.

**PRODUCT DESCRIPTION AND PACKAGING**  
**CitraFlow™** is supplied as a clear, plastic, prefilled syringe, packaged in a heat-sealed pouch. Each **CitraFlow™** prefilled syringe contains a sterile, clear and colorless solution of 0.467 g/mL sodium citrate; pH is adjusted with citric acid. **CitraFlow™** is latex free and contains no preservatives.

Catalogue Number	Description	Packaging
38143	3 mL Sodium Citrate Solution 46.7% (w/v) in a 5 mL syringe	One (1) plastic syringe packaged in a heat sealed pouch; 150 pouches (150 syringes) per product case
38143-1	3 mL Sodium Citrate Solution 46.7% (w/v) in a 5 mL syringe	Two (2) plastic syringes packaged in a heat sealed pouch; 100 pouches (200 syringes) per product case

**STORAGE AND STABILITY**  
**CitraFlow™** should be stored between 15°C and 30°C (59°F and 86°F) and protected from direct sunlight. **CitraFlow™** may be stored up to the expiration date on the product package.

**SYMBOLS ON PRODUCT LABELS**

REF: Catalogue (Reference) Number	Caution, consult accompanying documents	STERILE R Sterile by irradiation
LOT Lot number	Does not contain latex	Do not reuse, single use only
Use by	Recommended Storage Temperature	Authorised representative in the Emergo-Europe Prinsessegracht 20 2514 AP, The Hague The Netherlands
CE marking	Fill volume	EC REP

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