CitraFlow Plus™

Prefilled Syringe with Anticoagulant Sodium Citrate 4% (w/v) USP and Antimicrobial Ethanol 30% (v/v) Solution
Sterile – Non Pyrogenic – Single Use

INDIRECT IN-LINE ANTICOAGULANT AND ANTIMICROBIAL SOLUTION**INDICATIONS AND CLINICAL USE**

CitraFlow Plus™ is indicated for use in maintaining patency and preventing microbial growth in Vascular Access Devices (VAD's).

CONTRAINDICATIONS

CitraFlow Plus™ should not be used in patients with allergies or hypersensitivity to Sodium Citrate or Ethanol or who have liver disease, a history of alcohol abuse, weigh less than 20 kg or in pregnant women or nursing mothers. **CitraFlow Plus™** is not intended as a primary treatment for infected catheters.

WARNINGS AND PRECAUTIONS

Use aseptic technique. Single use only. Do not reuse. Not for direct intravenous injection. Do not use if package is opened or damaged. Do not use if the cap on the syringe is not intact. Visually inspect the contents of each prefilled syringe for clarity, particulate matter, precipitate, discoloration and leakage prior to use. Solutions showing haziness, particulate matter, precipitate, discoloration or leakage should not be used. Use with alcohol-compatible VAD.

ADVERSE REACTIONS

No known adverse reactions have been reported when the product is used as indicated.

INSTRUCTIONS FOR USE**USE ASEPTIC TECHNIQUE.**

1. Prior to initiation of a venous access session, aspirate and discard the solution that was instilled in the Vascular Access Device (VAD) during the previous session.
2. At the completion of each venous access session, flush the lumens of the VAD with sterile 0.9% Sodium Chloride USP solution to remove any blood, medication or other substance remaining in the VAD.
3. Aseptically open package containing **CitraFlow Plus™** prefilled syringe(s). Inspect the contents of the syringe(s) according to WARNINGS AND PRECAUTIONS.
4. Remove and discard syringe cap. Expel air from syringe. Do not allow air to be trapped in fluid path. Connect the syringe to the VAD. Use only with compatible luer lock connectors.
5. Slowly, over a 5 to 10 second period, instill the exact priming volume of **CitraFlow Plus™** solution into each VAD lumen, as a locking agent, in volumes corresponding to the luminal capacity as specified by the VAD manufacturer.
6. Remove and discard unused portions and empty syringes according to the institution's biohazardous waste policy.

ACTION

The **CitraFlow Plus™** solution acts as an indirect inline anticoagulant and antimicrobial agent for Vascular Access Device (VAD) locking. Sodium Citrate Solution exerts its anticoagulant effect via the chelation of ionized Calcium in the blood and tissues by the Citrate ion which prevents activation of Calcium dependent pro coagulants. In the extracorporeal circuit, Sodium Citrate causes anticoagulation by chelating ionized Calcium into a soluble complex. Calcium is an integral ion involved in the clotting cascade. Local removal by Citrate prevents the activation of clotting cofactors, factor X, and prothrombin and the ultimate formation of fibrin. Ethanol exerts its antimicrobial effect, via the denaturation of microorganisms, preventing biofilm formation which could result in VAD related blood stream infections. When used as indicated, **CitraFlow Plus™** eliminates the risks of systemic anticoagulation.

PRODUCT DESCRIPTION AND PACKAGING

CitraFlow Plus™ is supplied as a clear, plastic, prefilled syringe, packaged in a heat-sealed pouch. Each **CitraFlow Plus™** prefilled syringe contains a sterile, clear and colorless solution of 0.04 g/mL Sodium Citrate Dihydrate USP and 0.3 mL/mL of Ethanol; pH adjusted with citric acid. **CitraFlow Plus™** is latex free, non-pyrogenic, and contains no preservatives.

Catalogue Number	Description	Packaging
3844	5 mL Sodium Citrate 4% (w/v) USP with Ethanol 30% (v/v) Solution in a 5 mL syringe	One (1) plastic syringe packaged in a heat sealed pouch; 120 pouches (120 syringes) per product case
38443-1	3 mL Sodium Citrate 4% (w/v) USP with Ethanol 30% (v/v) Solution in a 5 mL syringe	Two (2) plastic syringes packaged in a heat sealed pouch; 100 pouches (200 syringes) per product case
38443	3 mL Sodium Citrate 4% (w/v) USP with Ethanol 30% (v/v) Solution in a 5 mL syringe	One (1) plastic syringe packaged in a heat sealed pouch; 150 pouches (150 syringes) per product case
38453	3 mL Sodium Citrate 4% (w/v) USP with Ethanol 30% (v/v) Solution in a 10 mL syringe	One (1) plastic syringe packaged in a heat sealed pouch; 100 pouches (100 syringes) per product case
38455	5 mL Sodium Citrate 4% (w/v) USP with Ethanol 30% (v/v) Solution in a 10 mL syringe	One (1) plastic syringe packaged in a heat sealed pouch; 100 pouches (100 syringes) per product case

STORAGE AND STABILITY

CitraFlow Plus™ should be stored between 15°C and 30°C (59°F and 86°F) and protected from direct sunlight. **CitraFlow Plus™** may be stored as indicated by the expiration date on the product package.

SYMBOLS ON PRODUCT LABELS

 Sterile by irradiation	REF: Catalogue (Reference) number	 Lot number
 Caution, consult accompanying documents	 Do not reuse, single use only	 Does not contain latex
 Use by	 Recommended Storage Temperature	 Fill volume (ml):

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