**PraxiFlow™**  
 1 Prefilled Syringe with Normal Saline Solution  
 1 Prefilled Syringe with Sodium Citrate Solution

**INDICATIONS AND CLINICAL USE**  
**Praxiject™** Normal Saline is for vascular access devices (VADs) flushing. **CitraFlow™** is indicated for use in maintaining patency of VADs.

**CONTRAINDICATIONS**  
 Patients with known allergies or hypersensitivity to constituents. For a listing see the PRODUCT DESCRIPTION AND PACKAGING section. Not for medication reconstitution or dilution.

**WARNINGS AND PRECAUTIONS**  
 Use aseptic technique. Single use only. Do not reuse. Not for direct intravenous injection. Do not use if package is opened or damaged. Do not use if the cap on the syringe is not intact. Visually inspect the contents of each prefilled syringe for clarity, particulate matter, precipitate, discoloration and leakage prior to use. Do not use if any of the aforementioned is observed. Do not allow air to be trapped in fluid path. Do not resterilize. Sterile unless package opened or damaged. May be placed on a sterile field.

**ADVERSE REACTIONS**  
 No known adverse reactions have been reported when the product is used as indicated. Parasthesia and/or dysgeusia may occur with high concentrations [30% & 46.7%] of Sodium Citrate for the product unintentionally passes into the vein.

**INSTRUCTIONS FOR USE**  
**USE ASEPTIC TECHNIQUE.**

- Prior to initiation of a venous access session, aspirate and discard, the solution that was instilled in the Vascular Access Device (VAD) during the previous session in accordance with institution's protocol.
- At the completion of each venous access session, aseptically open the package and inspect the contents of the syringes according to WARNINGS AND PRECAUTIONS.
- Remove and discard the **Praxiject™** Normal Saline syringe cap. Expel air from syringe. Do not allow air to be trapped in fluid path. Use only with compatible luer lock connectors. Flush the lumens of the VAD with the sterile saline solution to remove any blood, medication or other substance remaining in the VAD.
- Remove and discard the **CitraFlow™** syringe cap. Expel air from syringe. Do not allow air to be trapped in fluid path. Adjust to correct volume. Connect the syringe to the VAD. Use only with compatible luer lock connectors.
- Slowly, over a 5 to 10 second period, instill the exact priming volume of **CitraFlow™** solution into each VAD lumen, as a locking agent, in volumes corresponding to the luminal capacity as specified by the VAD manufacturer.
- Remove and discard unused portions and empty syringes according to the institution's biohazardous waste policy.

**ACTION**  
**Praxiject™** Normal Saline is a sterile, aqueous solution having approximately the same osmotic pressure and composition as extracellular fluids. It is non-irritating to tissues. It is used to flush VADs in order to maintain catheter patency and to prevent contact between incompatible medications or fluids.

**CitraFlow™** acts as an indirect in-line anticoagulant agent for VAD locking. In the extracorporeal circuit, Sodium Citrate causes anticoagulation by chelating ionized Calcium into a soluble complex. Calcium is an integral ion involved in the clotting cascade. Local removal by Citrate prevents the activation of clotting cofactors, factor X and prothrombin and the ultimate formation of fibrin. Systemic anticoagulation does not occur. The antimicrobial effect of **CitraFlow™** 30% and 46.7% is through binding and removal of Ca<sup>2+</sup> in the surrounding area. Ca<sup>2+</sup> may regulate certain genes responsible for growth and survival of microbes.

**PRODUCT DESCRIPTION AND PACKAGING**  
 Each kit is supplied as two (2) clear, plastic, prefilled syringes packaged in a heat sealed pouch. Each **Praxiject™** prefilled syringe contains a clear, colourless, sterile Normal Saline 0.9% solution. Each **CitraFlow™** prefilled syringe contains a clear, colourless, sterile solution of 0.04 g/ml sodium citrate and ethanol (excipient) [CitraFlow™ 4%], or 0.3 g/ml sodium citrate [CitraFlow™ 30%], or 0.467 g/ml sodium citrate [CitraFlow™ 46.7%]. pH is adjusted with citric acid. The kits are latex free and contain no preservatives.

Catalogue Number	Description	Packaging
3255	10 mL Praxiject™ 0.9% Normal Saline in a 10 mL syringe + 3 mL CitraFlow™ 4% (w/v) in a 5 mL syringe	One (1) 3705C and one (1) 3854J; 75 pouches per case
3225	10 mL Praxiject™ 0.9% Normal Saline in a 10 mL syringe + 3 mL CitraFlow™ 30% (w/v) in a 5 mL syringe	One (1) 3705C and one (1) 3824J; 75 pouches per case
3215	10 mL Praxiject™ 0.9% Normal Saline in a 10 mL syringe + 3 mL CitraFlow™ 46.7% (w/v) in a 5 mL syringe	One (1) 3705C and one (1) 3814J; 75 pouches per case

**STORAGE AND STABILITY**  
**PraxiFlow™** kits should be stored between 15°C and 30°C (59°F and 86°F) and protected from direct sunlight and freezing. The product may be stored up to the expiration date on the product package.

**SYMBOLS ON PRODUCT LABELS**

REF: Catalogue (Reference) Number	Caution, consult accompanying documents	Sterile by irradiation
Lot number	Does not contain latex	Do not reuse, single use only
Use by	Recommended Storage Temperature	Authorised representative in the European community Emergo-Europe Prinseesgracht 20 2514 AP, The Hague The Netherlands
CE 0413 CE marking	Fill volume	

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