

PraxiFill™

Automated syringe filler

“Fast, accurate and repeatable.
Syringe filling made easy!”

PraxiFill™ is an automated syringe filler intended to increase productivity of the syringe filling process by eliminating the tedious manual task of aspirating solution into the syringe and capping the filled syringe.

PraxiFill™ is designed to be used with MedXL pre-capped sterile syringes for increased efficiency and productivity.


3cc, 5cc, 10cc, 20cc and 50cc syringes are available.

Filling speed is up to 8 syringes/minute depending on syringe size.

Increases productivity

Saves time and money

Reduces repetitive motion injuries



Features & Benefits

- ✓ Fast, accurate, repeatable filling and capping.
- ✓ No need for calibration, auto zero on all functions.
- ✓ Small footprint (compact).
- ✓ Can fill 3cc, 5cc, 10cc, 20cc, and 50cc syringes.
- ✓ Simple and fast change over between different sizes.
- ✓ Minimum fill level from 10% of syringe nominal capacity.
- ✓ Filling accuracy and repeatability to +/- 1%.
- ✓ No drip filling feature.
- ✓ Caps are designed for freezing and thawing.



For more information,
please contact your
representative.