

PraxiFill™

Automated syringe filler


“ Fast, accurate and repeatable.
Syringe filling made easy! ”

PraxiFill™ is an automated syringe filler intended to increase productivity of the syringe filling process by eliminating the tedious manual task of aspirating solution into the syringe and capping the filled syringe.

PraxiFill™ is designed to be used with MedXL pre capped sterile syringes for increased efficiency and productivity.

3cc, 5cc, 10cc, 20cc and 50cc syringes are available.

Filling speed is up to 8 syringes/minute depending on syringe size.



Increases productivity

Saves time and money

Reduces repetitive motion injuries

Features & Benefits

- ✓ Fast, accurate, repeatable filling and capping.
- ✓ No need for calibration, auto zero on all functions.
- ✓ Small footprint (compact).
- ✓ Can fill 3cc, 5cc, 10cc, 20cc and 50cc syringes.
- ✓ Simple and fast change over between different sizes.
- ✓ Minimum fill level from 10% of syringe nominal capacity.
- ✓ Filling accuracy and repeatability to +/- 1%.
- ✓ No drip filling feature.
- ✓ Caps are designed for freezing and thawing.

Photo representative only. Product may not be exactly as shown.

Fs-VM-13/03e ©2019



For more information,
please contact your
representative.

 MedXL Inc.
Pointe-Claire (Qc) Canada H9R 1A3

 medxl.com
 +1.514.695.7474

